

# Anti-MonoMethyl-Histone H4 (Lys20) [R04-6U-3] Monoclonal Antibody

## AGMB05371

### Description

---

This Anti-MonoMethyl-Histone H4 (Lys20) [R04-6U-3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                     |
|-----------------------|-----------------------------------------------------|
| <b>SKU:</b>           | AGMB05371                                           |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml) |
| <b>Synonyms:</b>      | H4K20me, H4, H4/n, H4F2, H4FN, FO108, HIST2H4       |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-F</b> <b>IHC-P</b> <b>ICC/IF</b>   |
| <b>Research Area:</b> | Epigenetics                                         |
| <b>Form:</b>          | Liquid                                              |

### Antibody Data

---

|                      |                     |
|----------------------|---------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat   |
| <b>Clone:</b>        | R04-6U-3            |
| <b>Clonality:</b>    | Monoclonal Antibody |
| <b>SwissProt ID:</b> | P62805              |
| <b>Immunogen:</b>    | Peptide             |
| <b>Gene ID:</b>      | -                   |
| <b>Gene Name:</b>    | H4C1                |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                           |
|--------------------------|-------------------------------------------|
| <b>Host Species:</b>     | Rabbit                                    |
| <b>Isotype:</b>          | IgG                                       |
| <b>Purification:</b>     | Affinity Purified                         |
| <b>Conjugated:</b>       | Unconjugated                              |
| <b>Modification:</b>     | Monomethylated                            |
| <b>Molecular Weight:</b> | Calculated MW:11 kDa, Observed MW: 11 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b>           | Application                                                                                                                                                                                                                                             | Antibody Dilution Ratio |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|
|                                           | WB                                                                                                                                                                                                                                                      | 1:500-1:1000            |
|                                           | IHC                                                                                                                                                                                                                                                     | 1:50-1:100              |
|                                           | IF                                                                                                                                                                                                                                                      | 1:50-1:200              |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |                         |