

Anti-GAPDH [R03-5G5] Monoclonal Antibody

AGMB05386

Description

This Anti-GAPDH [R03-5G5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05386
Contents:	1mL, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Gapd, Peptidyl-cysteine S-nitrosylase GAPDH
Applications:	WB IHC-F IHC-P ICC/IF FC IP
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat, Monkey, Rabbit, Chicken, Zebrafish, Hamster, Escherichia Coli
Clone:	R03-5G5
Clonality:	Monoclonal Antibody
SwissProt ID:	P16858
Immunogen:	Recombinant protein of human GAPDH
Gene ID:	14433

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: GAPDH
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 36 kDa, Observed MW: 36 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:2000
	ICC/IF	1:200-1:500
	FC	1:100-1:200
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.