

Anti-Lamin B1 (9B10) [9B10-7C6-8A7] Monoclonal Antibody

AGMB05396

Description

This Anti-Lamin B1 (9B10) [9B10-7C6-8A7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05396
Contents:	1mL, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LMN, ADLD, LMN2, LMNB, lamin B1, Lamin-B1,
Applications:	WB IHC-P ICC/IF ELISA
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	9B10-7C6-8A7
Clonality:	Monoclonal Antibody
SwissProt ID:	P20700
Immunogen:	Recombinant Protein of Lamin-B1
Gene ID:	4001
Gene Name:	LMNB1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 66 kDa, Observed MW: 66 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:2000
	IHC	1:50-1:100
	IF	1:100-1:500
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	