

Anti-beta Tubulin (Zebrafish Specific) [6H2-4E6-8G2] Monoclonal Antibody

AGMB05421

Description

This Anti-beta Tubulin (Zebrafish Specific) [6H2-4E6-8G2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05421
Contents:	1mL, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Beta 4 tubulin, TBB5, TUBB2, TUBB2A, tubulin beta 2A, beta tubulin, tubulin beta chain, Tubulin beta 2b, Alpha tubulin, Alpha-tubulin ubiquitous,
Applications:	WB
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Zebrafish
Clone:	6H2-4E6-8G2
Clonality:	Monoclonal Antibody
SwissProt ID:	P07437
Immunogen:	Synthetic Peptide of β-Tubulin
Gene ID:	203068

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TUBB
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 50 kDa, Observed MW: 55 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:5000-1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.