

Anti-beta I Tubulin [8E6-1E5-10A8] Monoclonal Antibody

AGMB05422

Description

This Anti-beta I Tubulin [8E6-1E5-10A8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05422
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	beta I tubulin, TBB1, TUBB1, Tubulin beta 1, Class VI beta tubulin, 2810484G07Rik, beta 1 tubulin, Beta tubulin 1, class VI,
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Zebrafish
Clone:	8E6-1E5-10A8
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9H4B7
Immunogen:	Synthetic Peptide of β I tubulin
Gene ID:	81027

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TUBB1
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 50 kDa, Observed MW: 50 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.