

Anti-Cystatin C [7F6-A5-F3] Monoclonal Antibody, unconjugated, detector

AGMB05438

Description

This Anti-Cystatin C [7F6-A5-F3] Monoclonal Antibody, unconjugated, detector is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05438
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CST3, Cystatin-C, Cystatin-3, Gamma-trace, Neuroendocrine basic polypeptide, Post-gamma-globulin, AD 8, AD8, Amyloid angiopathy and cerebral hemorrhage,
Applications:	ELISA
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	7F6-A5-F3
Clonality:	Monoclonal Antibody
SwissProt ID:	P01034
Immunogen:	A synthesized peptide derived from CST3

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	1471
Gene Name:	CST3
Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.