

Anti-FABP3 [3E8-B10] Monoclonal Antibody, unconjugated, detector

AGMB05440

Description

This Anti-FABP3 [3E8-B10] Monoclonal Antibody, unconjugated, detector is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05440
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	422 protein, Cardiac Fatty Acid Binding Protein, FABP 11, FABP 3, FABP11, FABP3, FABPH_HUMAN, fatty acid binding protein 11, Fatty acid binding protein 3, Fatty acid binding protein 3 muscle and heart, Fatty acid binding protein 3 muscle and heart mammary derived growth inhibitor, Fatty acid binding protein 3 muscle, Fatty acid binding protein 3
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	3E8-B10
Clonality:	Monoclonal Antibody
SwissProt ID:	P05413
Immunogen:	A synthesized peptide derived from FABP3

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	2170
Gene Name:	FABP3
Host Species:	Mouse
Isotype:	IgG2b
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.