

Anti-IL-8 [R02-3Z9] Monoclonal Antibody,Capture (PBS Only)

AGMB05453

Description

This Anti-IL-8 [R02-3Z9] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05453
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	IL-8 (3A4), 9E3, CXCL8, CEF-4, GCP-1, Granulocyte chemotactic protein 1
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R02-3Z9
Clonality:	Monoclonal Antibody
SwissProt ID:	P10145
Immunogen:	Recombinant protein of human IL18
Gene ID:	3576
Gene Name:	CXCL8

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.