

Anti-Cytokeratin 19 [R07-5N1] Monoclonal Antibody,Capture (PBS Only) AGMB05455

Description

This Anti-Cytokeratin 19 [R07-5N1] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05455
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CK-19, KRT19, Keratin, type I cytoskeletal 19, Cytokeratin-19, Keratin-19, K19
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-5N1
Clonality:	Monoclonal Antibody
SwissProt ID:	P08727
Immunogen:	Recombinant protein of human Cytokeratin 19
Gene ID:	3880
Gene Name:	KRT19

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.