

Anti-CCL18 [R05-6B4] Monoclonal Antibody,Capture (PBS Only)

AGMB05457

Description

This Anti-CCL18 [R05-6B4] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05457
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CKb7, PARC, AMAC1, DCCK1, MIP-4, AMAC-1, DC-CK1, SCYA18
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-6B4
Clonality:	Monoclonal Antibody
SwissProt ID:	P55774
Immunogen:	Recombinant protein of human CCL18
Gene ID:	6362
Gene Name:	CCL18

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.