

**Anti-IL-2 Receptor alpha [R09-7P1] Monoclonal Antibody,
Unconjugated, Detector (PBS Only)**
AGMB05492

Description

This Anti-IL-2 Receptor alpha [R09-7P1] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05492
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	IL2RA, Interleukin-2 receptor subunit alpha, IL-2 receptor subunit alpha, IL-2-RA, IL-2R subunit alpha, IL2-RA, TAC antigen, p55, CD25
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-7P1
Clonality:	Monoclonal Antibody
SwissProt ID:	P01589
Immunogen:	Recombinant protein
Gene ID:	3559

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	IL2RA
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.