

Anti-IL-3 [R09-1S1] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05494

Description

This Anti-IL-3 [R09-1S1] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05494
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	IL3, IL-3, MCGF, Interleukin-3, Hematopoietic growth factor, Mast cell growth factor (Mast cell growth factor,MCGF), Multipotential colony-stimulating factor, P-cell-stimulating factor
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-1S1
Clonality:	Monoclonal Antibody
SwissProt ID:	P08700
Immunogen:	Recombinant protein

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	3562
Gene Name:	IL3
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.