

Anti-beta Tubulin [R79-7O-5] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05570

Description

This Anti-beta Tubulin [R79-7O-5] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05570
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Beta 4 tubulin, TBB5, TUBB2, TUBB2A, tubulin beta 2A, beta tubulin, tubulin beta chain, Tubulin beta 2b, Alpha tubulin, Alpha-tubulin ubiquitous,
Applications:	ELISA
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R79-7O-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P07437
Immunogen:	Recombinant protein
Gene ID:	203068

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	TUBB
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage: SSStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.