

Anti-CD147 [R89-2G-1] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05572

Description

This Anti-CD147 [R89-2G-1] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05572
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	BSG, Basigin, 5F7, Collagenase stimulatory factor, Extracellular matrix metalloproteinase inducer, EMMPRIN, Leukocyte activation antigen M6, OK blood group antigen, Tumor cell-derived collagenase stimulatory factor, TCSF, CD antigen CD147
Applications:	ELISA
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R89-2G-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P35613
Immunogen:	Recombinant protein

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	682
Gene Name:	BSG
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.