

Anti-E Cadherin [R59-2N-9] Monoclonal Antibody,Capture (PBS Only) AGMB05585

Description

This Anti-E Cadherin [R59-2N-9] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05585
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Cdh1,CD324,Cadherin-1,Epithelial cadherin (E-cadherin),Uvomorulin,E-Cad/CTF1,E-Cad/CTF2,E-Cad/CTF3, Arc 1, CADH1_HUMAN, Cadherin 1,
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R59-2N-9
Clonality:	Monoclonal Antibody
SwissProt ID:	P09803
Immunogen:	Recombinant protein
Gene ID:	12550

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: Cdh1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.