

Anti-MIF [R85-5I-4] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05606

Description

This Anti-MIF [R85-5I-4] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05606
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Mif,GIF, DER6, Glif, macrophage migration inhibitory factor, Delayed early response protein 6 (DER6), Glycosylation-inhibiting factor, L-dopachrome tautomerase
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R85-5I-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P34884
Immunogen:	Recombinant protein

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	17319
Gene Name:	Mif
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage: Store at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.