

Anti-MMP12 [R47-4V-1] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05612

Description

This Anti-MMP12 [R47-4V-1] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05612		
Contents:	100µg Bradford Reagent: 1 vial (2ml)		
Synonyms:	Mme, Mmel,Mmp12,Macrophage metalloproteinase-12 (MMP-12)	metalloelastase,MME,Matrix	
Applications:	ELISA		
Research Area:	Cardiovascular		
Form:	Liquid		

Antibody Data

Reactivity:	Mouse
Clone:	R47-4V-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P34960
Immunogen:	Recombinant protein
Gene ID:	17381

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	Mmp12
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.