

Anti-MMP3 [R67-5G-2] Monoclonal Antibody,Capture (PBS Only)

AGMB05625

Description

This Anti-MMP3 [R67-5G-2] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05625		
Contents:	100µg Bradford Reagent: 1 vial (2ml)		
Synonyms:	Mmp3,Stromelysin-1,SL-1,EMS-2,Matrix 3),Transin-1	metalloproteinase-3	(MMP-
Applications:	ELISA		
Research Area:	Cardiovascular		
Form:	Liquid		

Antibody Data

Reactivity:	Rat
Clone:	R67-5G-2
Clonality:	Monoclonal Antibody
SwissProt ID:	P03957
Immunogen:	Recombinant protein
Gene ID:	171045

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: Mmp3

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.