

Anti-MMP8 [R11-8H-7] Monoclonal Antibody,Capture (PBS Only)

AGMB05629

Description

This Anti-MMP8 [R11-8H-7] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05629
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Mmp8,Neutrophil collagenase,Collagenase 2,Matrix metalloproteinase-8 (MMP-8)
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R11-8H-7
Clonality:	Monoclonal Antibody
SwissProt ID:	O70138
Immunogen:	Recombinant protein
Gene ID:	17394

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: Mmp8

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.