

Anti-NG2 [R53-2Z-1] Monoclonal Antibody,Capture (PBS Only)

AGMB05645

Description

This Anti-NG2 [R53-2Z-1] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05645
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	An2, Kiaa4232, Ng2,Cspg4,Chondroitin sulfate proteoglycan 4,Chondroitin sulfate proteoglycan NG2,Proteoglycan AN2
Applications:	ELISA
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R53-2Z-1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8VHY0
Immunogen:	Recombinant protein
Gene ID:	121021

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	Cspg4
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage: SSStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.