

Anti-P75 Neurogranin [R35-2F-4] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05658

Description

This Anti-P75 Neurogranin [R35-2F-4] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05658
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	neurogranin (protein kinase C substrate, RC3), RC3, hng
Applications:	ELISA
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R35-2F-4
Clonality:	Monoclonal Antibody
SwissProt ID:	Q92686
Immunogen:	Recombinant protein
Gene ID:	4900

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NRGN
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.