

# Anti-PDGFR alpha [R35-9P-2] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05662

## Description

---

This Anti-PDGFR alpha [R35-9P-2] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB05662
<b>Contents:</b>	100µg Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	PDGFRA, PDGFR2, RHEPDGFRA, Platelet-derived growth factor receptor alpha, PDGF-R-alpha, PDGFR-alpha, Alpha platelet-derived growth factor receptor, Alpha-type platelet-derived growth factor receptor, CD140 antigen-like family member A, CD140a antigen, Platelet-derived growth factor alpha receptor, Platelet-derived growth factor receptor 2, PDGFR-2, CD140a
<b>Applications:</b>	<b>ELISA</b>
<b>Research Area:</b>	Cardiovascular
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R35-9P-2
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P16234

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

<b>Immunogen:</b>	Recombinant protein
<b>Gene ID:</b>	5156
<b>Gene Name:</b>	PDGFRA
<b>Host Species:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Purification:</b>	Affinity Purified
<b>Conjugated:</b>	Unconjugated
<b>Modification:</b>	Unmodified
<b>Molecular Weight:</b>	-

## Preparation & Storage

---

**Storage:** SStore at 4°C.Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 0.01M PBS, pH 7.2

**Antibody Dilution Ratio:** Optimize antibody dilution ratio according to assay requirements and sample type.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.