

Anti-UCH-L1/PGP9.5 [R37-3A-8] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05668

Description

This Anti-UCH-L1/PGP9.5 [R37-3A-8] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05668
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	UCHL1 / PGP9.5, UCHL1, B220, CD 45, CD45, cd45 antigen, ec3.1.3.48, GP 180, GP180, Human homolog of severe combined immunodeficiency due to PTPRC deficiency, L CA, L-CA, lca, Leukocyte common antigen, LY 5, LY5, Protein tyrosine phosphatase receptor type C, Protein tyrosine phosphatase receptor type c polypeptide, PTPRC, PTPRC_HUMAN, Receptor-type tyrosine-protein phosphatase C, SCID due to PTPRC deficiency, t200, T200 glycoprotein, T200 leukocyte common antigen.
Applications:	ELISA
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R37-3A-8
Clonality:	Monoclonal Antibody

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

SwissProt ID:	P09936
Immunogen:	Recombinant protein
Gene ID:	7345
Gene Name:	UCHL1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.