

Anti-S100B [R75-9M-9] Monoclonal Antibody,Capture (PBS Only)

AGMB05675

Description

This Anti-S100B [R75-9M-9] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05675
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	NEF, Protein S100 B, Protein S100-B, S 100 calcium binding protein beta chain, S 100 protein beta chain, S-100 protein beta chain, S-100 protein subunit beta, S100, S100 calcium binding protein beta (neural), S100 calcium-binding protein B, S100 protein beta chain, S100B, S100B_HUMAN, S100beta.
Applications:	ELISA
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R75-9M-9
Clonality:	Monoclonal Antibody
SwissProt ID:	P04271
Immunogen:	Recombinant protein

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 6285
Gene Name: S100B
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.