

Anti-S100A4 [R62-2M-9] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05682

Description

This Anti-S100A4 [R62-2M-9] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05682
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	18A2, 42A, calcium Placental protein, Calvasculin, CAPL, Fibroblast specific protein 1, Fibroblast specific protein, FSP1, Leukemia multidrug resistance associated protein, Malignant transformation suppression 1, Metastasin, MTS1, OTTHUMP00000015467, OTTHUMP00000015468, P9KA, PEL98, Placental calcium-binding protein, Protein Mts1, Protein S100 A4, Protein S100-A4, S100 calcium binding protein A4 (calcium protein, calvasculin, metastasin, murine placental homolog), S100 calcium binding protein A4, S100 calcium-binding protein A4, S100a4, S10A4_HUMAN.
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R62-2M-9
Clonality:	Monoclonal Antibody

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

SwissProt ID:	P26447
Immunogen:	Recombinant protein
Gene ID:	6275
Gene Name:	S100A4
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	Store at 4°C. Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.