

Anti-S100A9 [R19-5C-9] Monoclonal Antibody,Capture (PBS Only)

AGMB05685

Description

This Anti-S100A9 [R19-5C-9] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05685
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	S100 alpha9, S100A9, Calgranulin B, CFAG, CAGB, Calgranulin-B, CGLB, L1Ag, MRP14, NIF, Protein S100-A9, LIAG, MAC387, MRP-14
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R19-5C-9
Clonality:	Monoclonal Antibody
SwissProt ID:	P06702
Immunogen:	Recombinant protein
Gene ID:	6280

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: S100A9

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.