

Anti-VEGF Receptor 3 [R35-6N-6] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05714

Description

This Anti-VEGF Receptor 3 [R35-6N-6] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05714
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Flt-4, Vegfr3, Flt4, Vascular endothelial growth factor receptor 3, VEGFR-3, Fms-like tyrosine kinase 4 (FLT-4), Tyrosine-protein kinase receptor FLT4
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R35-6N-6
Clonality:	Monoclonal Antibody
SwissProt ID:	P35917
Immunogen:	Recombinant protein
Gene ID:	14257

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	Flt4
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.