

Anti-YKL-40/CHI3L1 [R34-4Z-7] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05722

Description

This Anti-YKL-40/CHI3L1 [R34-4Z-7] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05722
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Brp39, Chil1,Chi3l1,Chitinase-3-like protein 1,BRP39 protein,Breast regression protein 39,Cartilage glycoprotein 39 (CGP-39, GP-39), 39 kDa synovial protein, ASRT7, Cartilage glycoprotein 39,
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R34-4Z-7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q61362
Immunogen:	Recombinant protein

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	12654
Gene Name:	Chi3l1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	Store at 4°C. Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.