

**Anti-IL-6 [R09-1H-5] Monoclonal Antibody,unconjugated,
Unconjugated, Detector (PBS Only)
AGMB05738**

Description

This Anti-IL-6 [R09-1H-5] Monoclonal Antibody,unconjugated, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05738
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	IL6, IFNB2, Interleukin-6, IL-6, B-cell stimulatory factor 2, BSF-2, CTL differentiation factor, CDF, Hybridoma growth factor, Interferon beta-2, IFN-beta-2
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-1H-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P05231
Immunogen:	Recombinant protein

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	3569
Gene Name:	IL6
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage: Store at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.