

# Anti-AMBP [R03-2X-9] Monoclonal Antibody,Capture (PBS Only)

## AGMB05741

### Description

---

This Anti-AMBP [R03-2X-9] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                          |
|-----------------------|----------------------------------------------------------|
| <b>SKU:</b>           | AGMB05741                                                |
| <b>Contents:</b>      | 100µg<br>Bradford Reagent: 1 vial (2ml)                  |
| <b>Synonyms:</b>      | A1M, HCP, ITI, UTI, EDC1, HI30, ITIL, IATIL, ITILC, AMBP |
| <b>Applications:</b>  | <b>ELISA</b>                                             |
| <b>Research Area:</b> | Cardiovascular                                           |
| <b>Form:</b>          | Liquid                                                   |

### Antibody Data

---

|                      |                     |
|----------------------|---------------------|
| <b>Reactivity:</b>   | Human               |
| <b>Clone:</b>        | R03-2X-9            |
| <b>Clonality:</b>    | Monoclonal Antibody |
| <b>SwissProt ID:</b> | P02760              |
| <b>Immunogen:</b>    | Recombinant protein |
| <b>Gene ID:</b>      | 259                 |
| <b>Gene Name:</b>    | AMBP                |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                       |
|--------------------------|-----------------------|
| <b>Host Species:</b>     | Rabbit                |
| <b>Isotype:</b>          | IgG                   |
| <b>Purification:</b>     | Affinity Purification |
| <b>Conjugated:</b>       | Unconjugated          |
| <b>Modification:</b>     | Unmodified            |
| <b>Molecular Weight:</b> | -                     |

## Preparation & Storage

---

|                                           |                                                                                                                                                                                                                                                         |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Storage:</b>                           | SSStore at 4°C.Avoid freeze/thaw cycles.<br>Store Bradford Reagent at Room Temperature for 1 Year.                                                                                                                                                      |
| <b>Storage Buffer:</b>                    | Liquid in 0.01M PBS, pH 7.2                                                                                                                                                                                                                             |
| <b>Antibody Dilution Ratio:</b>           | Optimize antibody dilution ratio according to assay requirements and sample type.                                                                                                                                                                       |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |