

Anti-CD14 [R02-8J-6] Monoclonal Antibody,Capture (PBS Only)

AGMB05763

Description

This Anti-CD14 [R02-8J-6] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05763
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CD14, Monocyte differentiation antigen CD14, Myeloid cell-specific leucine-rich glycoprotein, CD antigen CD14
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R02-8J-6
Clonality:	Monoclonal Antibody
SwissProt ID:	P08571
Immunogen:	Recombinant protein
Gene ID:	929

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	CD14
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purification
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage: SSStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.