

# Anti-HLA-DRA [R09-5F-2] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05802

## Description

---

This Anti-HLA-DRA [R09-5F-2] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

|                       |                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB05802                                                                                         |
| <b>Contents:</b>      | 100µg<br>Bradford Reagent: 1 vial (2ml)                                                           |
| <b>Synonyms:</b>      | MLRW, HLA-DRA1, MHC class II antigen DRA, HLA class II histocompatibility antigen, DR alpha chain |
| <b>Applications:</b>  | <b>ELISA</b>                                                                                      |
| <b>Research Area:</b> | Immunology                                                                                        |
| <b>Form:</b>          | Liquid                                                                                            |

## Antibody Data

---

|                      |                     |
|----------------------|---------------------|
| <b>Reactivity:</b>   | Human               |
| <b>Clone:</b>        | R09-5F-2            |
| <b>Clonality:</b>    | Monoclonal Antibody |
| <b>SwissProt ID:</b> | P01903              |
| <b>Immunogen:</b>    | Recombinant protein |
| <b>Gene ID:</b>      | 3122                |

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                       |
|--------------------------|-----------------------|
| <b>Gene Name:</b>        | HLA-DRA               |
| <b>Host Species:</b>     | Rabbit                |
| <b>Isotype:</b>          | IgG                   |
| <b>Purification:</b>     | Affinity Purification |
| <b>Conjugated:</b>       | Unconjugated          |
| <b>Modification:</b>     | Unmodified            |
| <b>Molecular Weight:</b> | -                     |

## **Preparation & Storage**

---

**Storage:**               SSStore at 4°C.Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:**   Liquid in 0.01M PBS, pH 7.2

**Antibody Dilution Ratio:**   Optimize antibody dilution ratio according to assay requirements and sample type.

**Protein Quantification (Optional):**   To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.