

Anti-IGFBP7 [R05-6F-3] Monoclonal Antibody,Capture (PBS Only)

AGMB05811

Description

This Anti-IGFBP7 [R05-6F-3] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05811
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	IGFBP7, MAC25, PSF, Insulin-like growth factor-binding protein 7, IBP-7, IGF-binding protein 7, IGFBP-7, IGFBP-rP1, MAC25 protein, PGI2-stimulating factor, Prostacyclin-stimulating factor, Tumor-derived adhesion factor, TAF
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-6F-3
Clonality:	Monoclonal Antibody
SwissProt ID:	Q16270
Immunogen:	Recombinant protein
Gene ID:	3490

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: IGFBP7
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purification
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.