

# Anti-LOXL2 [R04-2H-6] Monoclonal Antibody,Capture (PBS Only)

## AGMB05835

### Description

---

This Anti-LOXL2 [R04-2H-6] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB05835
<b>Contents:</b>	100µg Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	LOR2, LOX L2, LOXL2, Lysyl oxidase homolog 2, Lysyl oxidase like 2, WS9 14, Lysyl oxidase-like protein 2, Lysyl oxidase-related protein 2, Lysyl oxidase-related protein WS9-14
<b>Applications:</b>	<b>ELISA</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R04-2H-6
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q9Y4K0
<b>Immunogen:</b>	Recombinant protein
<b>Gene ID:</b>	4017

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** LOXL2  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purification  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** SStore at 4°C.Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 0.01M PBS, pH 7.2

**Antibody Dilution Ratio:** Optimize antibody dilution ratio according to assay requirements and sample type.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.