

Anti-MUC1/CA15-3 [R01-7T-1] Monoclonal Antibody,Capture (PBS Only) AGMB05839

Description

This Anti-MUC1/CA15-3 [R01-7T-1] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05839
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CD227, MUC1, PUM, Mucin-1, MUC-1, Breast carcinoma-associated antigen DF3, Carcinoma-associated mucin, Episialin, H23AG, Krebs von den Lungen-6, KL-6, PEMT, Peanut-reactive urinary mucin, PUM, Polymorphic epithelial mucin, PEM, Tumor-associated ep
Applications:	ELISA
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R01-7T-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P15941
Immunogen:	Recombinant protein

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	4582
Gene Name:	MUC1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purification
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	Store at 4°C. Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.