

Anti-Mucin 16 [R03-6J-9] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05842

Description

This Anti-Mucin 16 [R03-6J-9] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05842
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	MUC16, CA125, Mucin-16, MUC-16, Ovarian cancer-related tumor marker CA125, CA-125, Ovarian carcinoma antigen CA125
Applications:	ELISA
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-6J-9
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8WXI7
Immunogen:	Recombinant protein
Gene ID:	94025

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	MUC16
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purification
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.