

Anti-Neutrophil Elastase [R06-4K-2] Monoclonal Antibody, Unconjugated, Detector (PBS Only) AGMB05844

Description

This Anti-Neutrophil Elastase [R06-4K-2] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05844
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	ELA2, ELANE, Elastase 2, Medullasin, Neutrophil elastase, PMN E, SCN1, Bone marrow serine protease,
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R06-4K-2
Clonality:	Monoclonal Antibody
SwissProt ID:	P08246
Immunogen:	Recombinant protein
Gene ID:	1991

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ELANE
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purification
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.