

Anti-Osteopontin [R01-2W-1] Monoclonal Antibody,Capture (PBS Only) AGMB05849

Description

This Anti-Osteopontin [R01-2W-1] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05849
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	SPP1, BNSP, OPN, PSEC0156, Osteopontin, Bone sialoprotein 1, Nephropontin, Secreted phosphoprotein 1, SPP-1, Urinary stone protein, Uropontin
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R01-2W-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P10451
Immunogen:	Recombinant protein
Gene ID:	6696

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	SPP1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purification
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.