

Anti-Collagen Type 1 [R05-1J-4] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05868

Description

This Anti-Collagen Type 1 [R05-1J-4] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05868
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Alpha-1 type I collagen, Alpha-2 type I collagen
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-1J-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P02452/P08123
Immunogen:	Recombinant protein
Gene ID:	1277;1278

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: COL1A1,COL1A2

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purification

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.