

Anti-SPARC [R06-6W-4] Monoclonal Antibody,Capture (PBS Only)

AGMB05877

Description

This Anti-SPARC [R06-6W-4] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05877
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	AA517111, SPARC, ON, SPARC, Basement-membrane protein 40, BM-40, Osteonectin, ON, Secreted protein acidic and rich in cysteine
Applications:	ELISA
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R06-6W-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P09486
Immunogen:	Recombinant protein
Gene ID:	6678

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SPARC
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purification
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.