

Anti-Tartrate Resistant Acid Phosphatase [R01-2W-8] Monoclonal Antibody,Capture (PBS Only)

AGMB05881

Description

This Anti-Tartrate Resistant Acid Phosphatase [R01-2W-8] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05881
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	HPAP, TRAP, TRAcP, TRACP5a, TRACP5b, TrATPase, ACP5
Applications:	ELISA
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R01-2W-8
Clonality:	Monoclonal Antibody
SwissProt ID:	P13686
Immunogen:	Recombinant protein
Gene ID:	54

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	ACP5
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purification
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage: SSStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.