

Anti-TIM3 [R04-9E-9] Monoclonal Antibody,Capture (PBS Only)

AGMB05917

Description

This Anti-TIM3 [R04-9E-9] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05917
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CD366, KIM-3, SPTCL, TIMD3, Tim-3, TIMD-3, HAVcr-2, TIM-3/HAVCR2
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R04-9E-9
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8TDQ0
Immunogen:	Recombinant protein
Gene ID:	84868
Gene Name:	HAVCR2

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purification
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.