

Anti-ICAM1 [4F5-F9-F13] Monoclonal Antibody, unconjugated, detector AGMB05940

Description

This Anti-ICAM1 [4F5-F9-F13] Monoclonal Antibody, unconjugated, detector is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05940
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CD antigen CD54, Intercellular adhesion molecule 1, ICAM-1, Major group rhinovirus receptor,
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	4F5-F9-F13
Clonality:	Monoclonal Antibody
SwissProt ID:	P05362
Immunogen:	Purified recombinant human CD54(ICAM-1) protein fragments expressed in E.coli.
Gene ID:	3383

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	ICAM1
Host Species:	Mouse
Isotype:	IgG2b
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.