

Anti-CGB3 [15B3] Monoclonal Antibody, unconjugated, detector

AGMB05944

Description

This Anti-CGB3 [15B3] Monoclonal Antibody, unconjugated, detector is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05944
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CGB, CGB3, CGB5, CGB7, CGB8, hCGB
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	15B3
Clonality:	Monoclonal Antibody
SwissProt ID:	P0DN86
Immunogen:	A synthesized peptide derived from native human chorionic gonadotropin beta
Gene ID:	94115,93659,1082

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	CGB3
Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in PBS containing 0.09% sodium azide as a preservative.
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.