

Anti-COL18A1/Endostatin [12G3] Monoclonal Antibody, Capture

AGMB05947

Description

This Anti-COL18A1/Endostatin [12G3] Monoclonal Antibody, Capture is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05947
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	KS, KNO, GLCC, KNO1, COL18A1, Antiangiogenic agent, COIA1_HUMAN, Col18a1,
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	12G3
Clonality:	Monoclonal Antibody
SwissProt ID:	P39060
Immunogen:	A synthesized peptide derived from endostatin
Gene ID:	80781

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	COL18A1
Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.