

# Anti-CEACAM5 [9H4] Monoclonal Antibody, unconjugated, detector

## AGMB05980

### Description

---

This Anti-CEACAM5 [9H4] Monoclonal Antibody, unconjugated, detector is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB05980
<b>Contents:</b>	100µg Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	CEA, CD66e, CEACAM5, Carcinoembryonic antigen-related cell adhesion molecule 5, Carcinoembryonic antigen, Meconium antigen 100
<b>Applications:</b>	<b>ELISA</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	9H4
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P06731
<b>Immunogen:</b>	A synthesized peptide derived from CEACAM5
<b>Gene ID:</b>	1048

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

<b>Gene Name:</b>	CEACAM5
<b>Host Species:</b>	Mouse
<b>Isotype:</b>	IgG
<b>Purification:</b>	Affinity Purified
<b>Conjugated:</b>	Unconjugated
<b>Modification:</b>	Unmodified
<b>Molecular Weight:</b>	-

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 0.09% sodium azide as a preservative.

**Antibody Dilution Ratio:** Optimize antibody dilution ratio according to assay requirements and sample type.

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.