

# Anti-ENO1 [9F10] Monoclonal Antibody, Capture

## AGMB05981

### Description

---

This Anti-ENO1 [9F10] Monoclonal Antibody, Capture is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                  |
|-----------------------|--------------------------------------------------|
| <b>SKU:</b>           | AGMB05981                                        |
| <b>Contents:</b>      | 100µg<br>Bradford Reagent: 1 vial (2ml)          |
| <b>Synonyms:</b>      | NNE, PPH, MPB1, ENO1L1, ENO1-IT1, HEL-S-17, ENO1 |
| <b>Applications:</b>  | <b>ELISA</b>                                     |
| <b>Research Area:</b> | Signal Transduction                              |
| <b>Form:</b>          | Liquid                                           |

### Antibody Data

---

|                      |                                         |
|----------------------|-----------------------------------------|
| <b>Reactivity:</b>   | Human                                   |
| <b>Clone:</b>        | 9F10                                    |
| <b>Clonality:</b>    | Monoclonal Antibody                     |
| <b>SwissProt ID:</b> | P06733                                  |
| <b>Immunogen:</b>    | A synthesized peptide derived from ENO1 |
| <b>Gene ID:</b>      | 2023                                    |
| <b>Gene Name:</b>    | ENO1                                    |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                   |
|--------------------------|-------------------|
| <b>Host Species:</b>     | Mouse             |
| <b>Isotype:</b>          | IgG1              |
| <b>Purification:</b>     | Affinity Purified |
| <b>Conjugated:</b>       | Unconjugated      |
| <b>Modification:</b>     | Unmodified        |
| <b>Molecular Weight:</b> | -                 |

## Preparation & Storage

---

|                                           |                                                                                                                                                                                                                                                         |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Storage:</b>                           | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.<br>Store Bradford Reagent at Room Temperature for 1 Year.                                                                                                      |
| <b>Storage Buffer:</b>                    | Liquid in PBS containing 0.09% sodium azide as a preservative.                                                                                                                                                                                          |
| <b>Antibody Dilution Ratio:</b>           | Optimize antibody dilution ratio according to assay requirements and sample type.                                                                                                                                                                       |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |