

Anti-KLK3 [3G9] Monoclonal Antibody, unconjugated, detector

AGMB05984

Description

This Anti-KLK3 [3G9] Monoclonal Antibody, unconjugated, detector is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05984
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	KLK3, APS, Prostate-specific antigen, PSA, Gamma-seminoprotein, Semin, Kallikrein-3, P-30 antigen, Semenogelase
Applications:	ELISA
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	3G9
Clonality:	Monoclonal Antibody
SwissProt ID:	P07288
Immunogen:	A synthesized peptide derived from KLK3
Gene ID:	354

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: KLK3
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.