

Anti-IgG2 [22C5] Monoclonal Antibody, unconjugated, detector

AGMB06018

Description

This Anti-IgG2 [22C5] Monoclonal Antibody, unconjugated, detector is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06018
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	IGHG2, Ig gamma-2 chain C region
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	22C5
Clonality:	Monoclonal Antibody
SwissProt ID:	P01859
Immunogen:	Recombinant human IgG2
Gene ID:	3501
Gene Name:	IGHG2

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in PBS containing 0.09% sodium azide as a preservative.
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.