

Anti-TIMP1 Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB06040

Description

This Anti-TIMP1 Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06040
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Metalloproteinase inhibitor 1, Erythroid-potentiating activity (EPA), Fibroblast collagenase inhibitor (Collagenase inhibitor), Tissue inhibitor of metalloproteinases 1 (TIMP-1), TIMP1, CLGI, TIMP
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	-
Clonality:	Monoclonal Antibody
SwissProt ID:	P01033
Immunogen:	-
Gene ID:	7076

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	TIMP1
Host Species:	Rabbit
Isotype:	-
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.